
基生研ゲノムインフォマティックス・トレーニングコース 2018夏 RNA-seq入門 - NGSの基礎からde novo解析まで-実践編: RNA-seq解析パイプライン 2018.07.26-2018.07.27

NGS基本ツールIGV

基礎生物学研究所 生物機能解析センター 山口勝司

データ可視化ツール・IGVの紹介・実習



なぜIGVを取り上げるか


データ可視化ツール

- 自分のパソコン(ローカル環境)にインストールして使うタイプ
- サーバーに構築して、ネットワーク上で使うタイプ

The JBrowse Genome Browser

JBrowse is a fast, embeddable genome browser built completely with JavaScript and HTML5, with optional run-once data formatting tools written in Perl.

Latest Release - JBrowse 1.14.2



後者はコミュニティーで広く利用、 あるいはウェブ公開を目的とするには 良いが、ネットワーク・情報セキュリティの 高度な知識も要求される。

より大容量なデータに対応できる。

管理者的な人がいて、その人がやって くれるなら、これも良いが・・・

もっとお手軽なものとしてIGVを紹介

可視化ツールに求められるものは何か

膨大なデータを如何に直感的に理解できるようにするか sortや絞り込みができる表データと対比双璧

多様なデジタル情報

- ·配列、GC ratio、遺伝子情報
- •遺伝子発現情報
- ·SNPの位置情報·頻度情報
- 様々なデータの精度情報
- ChIP-seq, RAD-seq, BS-seq

レファレンス配列 / gene model / gene annotationとNGSデータを並べて比較複数のデータセットを並べて比較

様々なスケールで比較・統合的に解釈できるようにしたい

ゲノムviewerに自分のデータを乗せ、 統合的直感的に比較・解釈できること

可視化ツールをどう選ぶか

選択の基準

genome data viewing に求められるもの 取捨選択の基準


- 1. 無料 / 有料 / 基本無料
- 2. 個人的レベルの使用 / コミュニティーレベルの使用
- 3. 見るだけ/自分から色々工夫
- 4.アクセスのしやすさ・使いやすさ 導入に必要なコンピュータスペック マニュアルは分かりやすいか 情報の多さ 利用の簡便さ 使っている人が近くにいるか

Integrative Genomics Viewer(IGV)


お手軽ツール


- アカデミックウェアで無料
- ・コミュニティーでの利用者が多いから、情報も多い
- ・javaのプログラムなので、オールプラットフォーム対応
- マニュアルは親切、サンプルデータのある
- •WEBサーバーではなく、PCレベルでできる
- ・データ閲覧環境の共有が可能


誰もが簡便に使えるものが良い。












Nature Biotech. 29:24-26 (2011) Supplement figureからの抜粋



Downloads

Integrative Genomics Viewer - IGV 2.4


Install IGV

NOTE: IGV 2.4.x requires Java 8. Java versions 9 and above are currently not supported.

Use one of the following 4 options to install and run the current version of IGV.

- Download and unzip the Mac App Archive, then double-click the IGV application to run it. The application can be
 moved to the Applications folder, or anywhere else.
- 2. Download and unzip the **Windows Zip Archive**, then double-click the *igv.bat* file to start IGV. A black console window will appear, followed by the IGV application. **Note**: Windows users with **high resolution screens** should use this version it includes a modified Java executable for use with high-resolution screens.
- 3. Download and unzip the Binary Distribution archive. IGV is launched from a command prompt follow the instructions in the *readme* file. To launch IGV on Mac or Linux use the shell script *igv.sh*. On Windows use *igv.bat*.
- 4. Click on one of the Launch buttons below to download a .jnlp file and execute the file using Java Web Start (JWS).
 - Mac users: If you are notified of security errors that prevent launching IGV, try the following:
 - Right-click on the downloaded .jnlp file; select Open With > Java Web Start; dismiss the warnings.
 - After IGV has been run this way at least once from the .jnlp file, you can double-click on the file to launch.
 - Windows users: To run with more than 1.2 GB of memory on Windows you must install 64-bit Java. Most Windows installs do not include 64-bit Java by default, even if the operating system is 64-bit. Attempting to use the 2GB or greater launch options with 32-bit Java will result in the error "could not create virtual machine".







ゲノムViewerなので次世代DNAシーケンサーのデータに限定されない。 マイクロアレイの結果や、ゲノムアノテーションの情報も随時表示できる。

対応するファイル形式に応じて、表示方法が決まる。




File Formats

IGV supports a number of different file formats for experimental data and genome annotations. For a complete list of supported formats see http://www.broadinstitute.org/igv/FileFormats. The following table shows the recommended file formats for a number of common data types.

Source Data	Recommended File Formats
ChIP-Seq, RNA-Seq	WIG, TDF
Copy number	CN, SNP, TDF, canary_calls (Birdsuite)
Gene expression data	GCT, RES, TDF
Genome annotations	GFF, BED, GTF, PSL, UCSC table format
GISTIC data	GISTIC
LOH data	LOH, TDF
Mutation data	MUT, MAF
Variant calls	VCF
RNAi data	GCT
Segmented data	SEG, CBS
Sequence alignment data	BAM, SAM, PSL
Any numeric data	IGV, WIG, TDF
Sample metatadata	Tab-delimited sample info file

公開情報のviewerとして



その他の便利機能

セッションの保存

表示しているデータの読み込み状況を、それごと保存。 セッションをロードすることで、意図した画面を表示できる。 データセットが揃っていること、フォルダー構造が同一である必要がある。


バッチ処理

重要領域の画面スナップショットを自動で取ったりできる。


new load myfile.bam snapshotDirectory mySnapshotDirectory genome hg18 goto chr1:65,289,335-65,309,335 sort position collapse snapshot goto chr1:113,144,120-113,164,120 sort base collapse snapshot

Inversions


An inversion is a large section of DNA that is reversed in the subject genome compared to the reference genome.



When an inversion shows up in paired-end reads, the reads are distinctively variant from the reference genome.



This appears in IGV as shown below



Interpreting Color by Insert Size

The inferred insert size can be used to detect structural variants, such as:


- deletions
- insertions
- inter-chromosomal rearrangements

IGV uses color coding to flag anomalous insert sizes. When you select Color alignments>by insert size in the popup menu, the default coloring scheme is:


- for an insert that is larger than expected
- for an insert that is smaller than expected
- 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 X Y for paired end reads that are coded by the chromosome on which their mates can be found

Deletions

A deletion is a large section of DNA that is absent in the subject genome compared to the reference genome.




The "expected" insert size is the insert size obtained in sequencing the subject genome. The "inferred" insert size is the insert size that would result in the reference genome, assuming the same pair of reads.




RNA-Seqのデータ表示させる





GCTファイルでgene ローカスの発現情報を図示





Gene listを定義して サンプルごと 条件ごと の発現・発現変動を カラーマップできる





IGV実習

Downloads Integrative Genomics Viewer - IGV 2.4 NOTE: IGV 2.4.x requires <u>Java 8</u>. Java versions 9 and above are currently not supported. Use one of the following 4 options to install and run the current version of IGV.

1. Download and unzip the Mac App Archive, then double-click the IGV application to run it. The application can be moved to the Applications folder, or anywhere else.

2. Download and unzip the **Windows Zip Archive**, then double-click the *igv.bat* file to start IGV. A black console window will appear, followed by the IGV application. **Note**: Windows users with **high resolution screens** should use this version — it includes a modified Java executable for use with high-resolution screens.

3. Download and unzip the Binary Distribution archive. IGV is launched from a command prompt – follow the instructions in the readme file. To launch IGV on Mac or Linux use the shell script igv.sh. On Windows use igv.bat.

4. Click on one of the Launch buttons below to download a .jnlp file and execute the file using Java Web Start (JWS).

- Mac users: If you are notified of security errors that prevent launching IGV, try the following:
 Right-click on the downloaded .jnlp file; select Open With > Java Web Start; dismiss the warnings.
 After IGV has been run this way at least once from the .jnlp file, you can double-click on the file to launch.
- Windows users: To run with more than 1.2 GB of memory on Windows you must install 64-bit Java. Most Windows installs do not include 64-bit Java by default, even if the operating system is 64-bit. Attempting to use the 2GB or greater launch options with 32-bit Java will result in the error "could not create virtual machine".

占 Launch

Launch with 750 MB

Launch with **1.2** GB (Max usable memory for Windows with 32-bit Java)

Launch with 2 GB (Max usable memory for 32-bit MacOS) Launch with 10 GB

(Only for large memory machines with 64-bit Java)

IGVの使用法を学ぶと共に 先のファイルフォーマットも 確認しよう

以下のファイルを確認


buc.genome.fasta


buc.gtf


buc_cg.wig


illumina_ex_B2_Read_bowtie2.mate.sort.bam illumina_ex_B2_Read_bowtie2.mate.sort.bam.bai illumina ex_B4_Read_bowtie2.mate.sort.bam illumina ex B4 Read bowtie2.mate.sort.bam.bai


31











































IGV紹介のまとめ

可視化ツールとして十分な機能を持つ

- -無料
- ・比較的簡単・お手軽
- ・自分で見るためにも良し、人に見せるためにも良し
- ・利用範囲は次世代DNAシーケンサーに限定しない 広くゲノミクスの解析に有用

ごく一部のみの機能を紹介しました。 ウェブサイトを見ながら復習をお勧めします。